

SHORT GENOME REPORT

Open Access



# Complete genome sequence of esterase-producing bacterium *Croceicoccus marinus* E4A9<sup>T</sup>

Yue-Hong Wu, Hong Cheng, Ying-Yi Huo, Lin Xu, Qian Liu, Chun-Sheng Wang and Xue-Wei Xu\*

## Abstract

*Croceicoccus marinus* E4A9<sup>T</sup> was isolated from deep-sea sediment collected from the East Pacific polymetallic nodule area. The strain is able to produce esterase, which is widely used in the food, perfume, cosmetic, chemical, agricultural and pharmaceutical industries. Here we describe the characteristics of strain E4A9, including the genome sequence and annotation, presence of esterases, and metabolic pathways of the organism. The genome of strain E4A9<sup>T</sup> comprises 4,109,188 bp, with one chromosome (3,001,363 bp) and two large circular plasmids (761,621 bp and 346,204 bp, respectively). Complete genome contains 3653 coding sequences, 48 tRNAs, two operons of 16S–23S–5S rRNA gene and three ncRNAs. Strain E4A9<sup>T</sup> encodes 10 genes related to esterase, and three of the esterases (E3, E6 and E10) was successfully cloned and expressed in *Escherichia coli* Rosetta in a soluble form, revealing its potential application in biotechnological industry. Moreover, the genome provides clues of metabolic pathways of strain E4A9<sup>T</sup>, reflecting its adaptations to the ambient environment. The genome sequence of *C. marinus* E4A9<sup>T</sup> now provides the fundamental information for future studies.

**Keywords:** *Croceicoccus marinus* E4A9<sup>T</sup>, Genome sequence, Esterase, Alphaproteobacteria

## Introduction

Lipolytic enzymes, including esterase (EC 3.1.1.1) and lipase (EC 3.1.1.3), are a general class of carboxylic ester hydrolases (EC 3.1.1), which catalyze the hydrolytic cleavage and formation of ester bonds [1, 2]. Esterase shows a preference for water-soluble short chain fatty acids (< 10 carbon atoms), while lipase prefers water-insoluble longer chain fatty acids (> 10 carbon atoms) [3, 4]. Many esterases do not require cofactors and have high stereospecificity toward chemicals, broad substrate specificity and high stability in organic solvents [4]. They are extensively used in the food, perfume, cosmetic, chemical, agricultural and pharmaceutical industries [5].

*Croceicoccus* [6], as a genus of the family *Erythrobacteraceae* [7], can be found in the marine environments, including deep-sea sediment, surface seawater and marine biofilm from a boat shell [6, 8, 9]. *C. marinus* E4A9<sup>T</sup>, the type strain of the genus *Croceicoccus*, was

isolated from deep-sea sediment collected from the East Pacific polymetallic nodule area [6]. The strain was able to produce esterase as well as lipase [6]. To get insight into the capability of esterase production, recently, we obtained the complete genome of *C. marinus* E4A9<sup>T</sup> and detected genes of esterase. This is the first genome report for the strain in the genus of *Croceicoccus*. We also describe the genomic sequencing related to its annotation for understanding their metabolic and ecological functions in the environment.


## Organism information

### Classification and features

*C. marinus* E4A9<sup>T</sup> was isolated from a deep-sea sediment sample collected from the East Pacific polymetallic nodule area (8°22'38" N, 145°23'56" W) at a depth of 5280 m (temperature 2 °C, salinity 3.4%). Strain E4A9<sup>T</sup> was obtained and routinely cultured on marine broth 2216 (MB, BD) at 30 °C. Subsequently polyphasic study of strain E4A9<sup>T</sup> was performed. A new species *Croceicoccus marinus* gen. nov. sp. nov. was proposed. Strain E4A9<sup>T</sup> is the type strain of the species of *C.*

\* Correspondence: xuxw@sio.org.cn

Key Laboratory of Marine Ecosystem and Biogeochemistry, Second Institute of Oceanography, State Oceanic Administration, 36th North BaoChu Road, Hangzhou 310012, China



**Fig. 1** Transmission electron microscopy showing the cell morphology (a) and ultrastructure (b) of *Croceicoccus marinus* E4A9<sup>T</sup>. The flagella are present. Bars represent scales of 0.5 μm (a) and 0.2 μm (b), respectively

*marinus* [6], and was deposited into the China General Microbiological Culture Collection (CGMCC 1.6776<sup>T</sup>).

*C. marinus* [6] is a valid species belonging to the family *Erythrobacteraceae* [7], in the order *Sphingomonadales* [10, 11], class *Alphaproteobacteria* [11, 12] and phylum *Proteobacteria* [13]. *C. marinus* E4A9<sup>T</sup> is a Gram-staining-negative and cocci-shaped bacterium (Fig. 1). It grew aerobically and used a series of organic carbon, such as L-arabinose, D-cellobiose, D-galactose and xylose, as sole sources of carbon and energy [6, 8]. Based on phylogenetic analysis of 16S rRNA gene sequence, the strain falls into the cluster comprising the *Croceicoccus* species with a high bootstrap value (Fig. 2). Interestingly, strain E4A9<sup>T</sup> could hydrolyze Tween 20, Tween 80 and tributyrin, indicating the presence of esterase as well as lipase [6]. The API ZYM system also supported the results that esterase (C4) and esterase lipase (C8) activities are present.

The general features of strain E4A9<sup>T</sup> was summarized in Table 1.


## Genome sequencing information

### Genome project history

*C. marinus* E4A9<sup>T</sup> [6] was selected for sequencing because it is relevant to genomic sequencing of the whole family of *Erythrobacteraceae* [7] and esterase production. The complete genome sequence was finished on May 29, 2015. The gap closure and annotation processes were performed by the authors. The GenBank accession number of the genome is CP019602, CP019603 and CP019604. The main genome sequence information is present in Table 2 and Table 3.

### Growth conditions and genomic DNA preparation

*C. marinus* E4A9<sup>T</sup> was aerobically cultivated in Marine Broth (MB, BD Difco™) at 30 °C and stored at -80 °C



**Fig. 2** Phylogenetic tree based on 16S rRNA gene sequences was constructed by neighbor-joining algorithms. Related sequences were aligned with Clustal W. Evolutionary distances were calculated according to the algorithm of the Kimura two-parameter model. Bootstrap values (> 60%) based on 1000 replications are shown at branch nodes. Filled circles indicate that the corresponding nodes were also recovered in the trees generated with the maximum-likelihood and maximum-parsimony algorithms. Bar, 0.01 substitutions per nucleotide position

**Table 1** Classification and general features of *Croceicoccus marinus* E4A9<sup>T</sup> according to the MIGS recommendations [30]

MIGS ID	Property	Term	Evidence code <sup>a</sup>
	Classification	Domain <i>Bacteria</i>	TAS [31]
		Phylum <i>Proteobacteria</i>	TAS [12]
		Class <i>Alphaproteobacteria</i>	TAS [11]
		Order <i>Sphingomonadales</i>	TAS [10]
		Family <i>Erythrobacteraceae</i>	TAS [7]
		Genus <i>Croceicoccus</i>	TAS [6]
		Species <i>Croceicoccus marinus</i> (Type) strain: Strain E4A9 <sup>T</sup> (CGMCC 1.6776 <sup>T</sup> = JCM 14846 <sup>T</sup> )	TAS [6]
		Gram stain	Negative
	Cell shape	Coccus	TAS [6]
	Motility	Motile	TAS [6]
	Sporulation	Non-sporulation	TAS [6]
	Temperature range	4–42 °C	TAS [6]
	Optimum temperature	28–30 °C	TAS [6]
	pH range; Optimum	6.0–9.0; 7.0	TAS [6]
	Carbon source	Organic carbon	TAS [6]
MIGS-6	Habitat	Deep-sea sediment	TAS [6]
MIGS-6.3	Salinity	Moderately halophilic, 0.5–10% NaCl	TAS [6]
MIGS-22	Oxygen requirement	Aerobic	TAS [6]
MIGS-15	Biotic relationship	Free-living	TAS [6]
MIGS-14	Pathogenicity	Non-pathogen	NAS
MIGS-4	Geographic location	East Pacific polymetallic nodule area	TAS [6]
MIGS-5	Sample collection	Not reported	
MIGS-4.1	Latitude	8°22'38" N	TAS [6]
MIGS-4.2	Longitude	145°23'56" W	TAS [6]
MIGS-4.4	Altitude	–5280 m	TAS [6]

<sup>a</sup>Evidence codes - *IDA* Inferred from Direct Assay, *TAS* Traceable Author Statement (i.e., a direct report exists in the literature), *NAS* Non-traceable Author Statement (i.e., not directly observed for the living, isolated sample, but based on a generally accepted property for the species, or anecdotal evidence). These evidence codes are from the Gene Ontology project [32]

**Table 2** Genome sequencing project information

MIGS ID	Property	Term
MIGS 31	Finishing quality	Finished
MIGS-28	Libraries used	10 kb
MIGS 29	Sequencing platforms	A PacBio RS II platform
MIGS 31.2	Fold coverage	248-fold
MIGS 30	Assemblers	HGAP Assembly version 2, Pacific Biosciences
MIGS 32	Gene calling method	GeneMarkS+ (NCBI)
	Locus Tag	A9D14
	Genbank ID	CP019602, CP019603, and CP019604
	GenBank Date of Release	June 13, 2017
	GOLD ID	Go0030822
	BIOPROJECT	PRJNA322659
MIGS 13	Source Material Identifier	CGMCC(China General Microbiological Culture Collection)
	Project relevance	Esterases production

**Table 3** Summary of genome: one chromosome and two plasmids

Label	Size (Mb)	Topology	INSDC identifier	RefSeq ID
Chromosome	3.001363	Linear	CP019602.1	NZ_CP019602.1
Plasmid 1 (pCME4A9I)	0.761621	Linear	CP019603.1	NZ_CP019603.1
Plasmid 2 (pCME4A9II)	0.346204	Linear	CP019604.1	NZ_CP019604.1

with 30% (*v/v*) glycerol. High-quality genomic DNA was extracted using the Qiagen DNA extraction kit, according to its protocol.

### Genome sequencing and assembly

The genome of strain E4A9<sup>T</sup> was sequenced using SMRT technology with a PacBio RS II platform (Zhejiang Tianke Co. Ltd., China). One library was constructed with 10 kb insert size according to the large SMRTbell gDNA protocol (Pacific Biosciences, USA). The sequencing generated 85,372 reads with an average length of 11,938 nt (972 Mb, 248-fold genome coverage). The de novo assembly of the reads was performed using HGAP Assembly version 2 (Pacific Biosciences, USA). The circularization of final contigs was checked and the overlapping ends were trimmed.

### Genome annotation

The rRNA genes were found via RNAmmer 1.2 Server [14] and tRNA genes were identified using tRNAscan-SE 2.0 online server [15]. The open reading frames (ORFs) and the functional annotation of translated ORFs were performed using the RAST server online [16] and GeneMarkS+. Classification of some predicted genes were analyzed using COG database [17] and Pfam [18]. Genes with signal peptides were predicted using SignalP 4.1 Server [19]. Genes with transmembrane helices were performed using TMHMM Server v. 2.0 [20]. The clustered regularly interspaced short palindromic repeats structures of the genomes were searched by CRISPRfinder program online [21]. Translated genes were assigned to Kyoto Encyclopedia of Genes and Genomes pathway using KEGG automatic annotation server with BBH method [22, 23]. The circular map of chromosome and plasmids were obtained using a CG View online server [24].

### Genome properties

The general features of strain E4A9 information are displayed in Table 1 and Table 2. The complete genome comprises 4,109,188 bp, with one chromosome (3,001,363 bp) and two large circular plasmids (plasmid pCME4A9I, 761,621 bp and plasmid pCME4A9II, 346,204 bp, respectively) (Fig. 3). The G + C content was 64.5 mol%. The genome of strain E4A9 contains 3653 coding sequences (CDSs), 48 tRNAs, two operons of 16S–23S–5S rRNA gene and three ncRNAs. Among the


genes, 132 were assigned to pseudogene. The summary of features and statistics of the genome is shown in Table 4 and genes belonging to COG functional categories are listed in Table 5.

Three replicons of the genome of strain E4A9, located in a circular chromosome and two large plasmids, were detected. Two plasmid replication initiator protein genes (ARU17925 and ARU18299) were found in the two plasmid sequence respectively, indicating that the genome of strain E4A9 contains two large circular plasmids. The G + C content of the two plasmids (63.5 mol% and 60.7 mol%, respectively) was a little lower than that of the chromosome (65.2 mol%). The two plasmids have high gene density with 702 and 303 protein-coding regions, respectively. Many unsuspected genes involved in metabolism of aromatic compounds were identified in plasmid pCME4A9I. Almost 10% of the plasmid pCME4A9II sequence carries genes encoding gene of subsystem feature virulence, disease and defense, and most of them were of the copper homeostasis and cobalt-zinc-cadmium resistance. The functions of these genes are consistent with the notion that the two plasmids play an important role in the adaptation of the bacteria in the sediment environment.

### Insights from the genome sequence

#### Esterases presence of *C. marinus* E4A9<sup>T</sup>

The presence of genes for the biotechnologically important enzymes like lipolytic enzymes were also predicted. Ten novel esterases were predicted (Fig. 4), and their amino acid sequences shared 58% to 85% identities to those of other lipolytic enzymes in the database. Phylogenetic analysis showed that predicted esterases E3 and E6 were grouped into family VII lipolytic enzymes and E10 was grouped into family II lipolytic enzymes. In order to investigate the biochemical properties of the esterases (E3, E6 and E10), recombinant plasmids were constructed and expressed in *Escherichia Coli* [25, 26]. After incubation of recombinant colonies for 48 h on the plate (Luria-Bertani agar medium) supplemented with 1% tributyrin, the three recombinant colonies had clear zones around the colonies. It indicated the presence of lipolytic activity. The calculated molecular weight of E3, E6 and E10 was 55.9, 46.1 and 22.4 kDa, respectively. The recombinant protein was soluble and purified using a Ni-NTA affinity chromatography column. The activity of purified E3, E6 and E10 was



**Fig. 3** Circular map of the chromosome (a), plasmid pCME4A9I (b) and plasmid pCME4A9II (c). From outside to the center: CDSs and RNA genes on the forward strand (colored by COG categories), CDSs and RNA genes on the reverse strand (colored by COG categories), G + C content (peaks out/inside the circle indicate values higher or lower than the average G + C content, respectively), GC skew (calculated as (G-C)/(G + C), green/purple peaks out/inside the circle indicate values higher or lower than 1, respectively), genome size

examined using *p*-nitrophenyl butyrate as substrate, and they had specific activities under standard reaction conditions (data not shown).

**Metabolism of *C. marinus* E4A9<sup>T</sup>**

The complete genome of *C. marinus* E4A9<sup>T</sup> was annotated for understanding the metabolic potentials based

on the key genes of metabolic pathways of carbon, nitrogen, sulfur and phosphorus. (i) Carbon metabolism. The genome of strain E4A9<sup>T</sup> is lack of carbon fixation and CO-oxidizing (*cox*) genes, indicating that the strain is not able to grow autotrophically. Strain E4A9<sup>T</sup> can use organic carbon sources (Table 1). The genome has a complete glycolysis pathway (Embden-Meyerhoff-Parnas

**Table 4** Genome statistics of *Croceicoccus marinus* E4A9<sup>T</sup>


Attribute	Value	% of Total
Genome size (bp)	4,109,188	100
DNA coding (bp)	3,565,753	86.78
DNA G + C (bp)	2,650,881	64.51
DNA scaffolds	3	–
Total genes	3842	100
Protein coding genes	3653	95.08
RNA genes	57	1.48
Pseudo genes	132	3.47
Genes in internal clusters	517	13.46
Genes with function prediction	2699	70.25
Genes assigned to COGs	2827	73.58
Genes with Pfam domains	1566	40.76
Genes with signal peptides	304	7.91
Genes with transmembrane helices	755	19.65
CRISPR repeats	1	0.03

pathway). In addition, it possesses key genes of the Entner-Doudoroff pathway, the pentose phosphate pathway, and the tricarboxylic acid cycle. (ii) Nitrogen metabolism. The genome of *C. marinus* E4A9<sup>T</sup> possesses ammonium transporter genes and amino acids transporter genes (e.g. methionine and L-proline/glycine betaine). Genes encoding enzymes involved in polyamines biosynthesis are present, but the lack of polyamines transporters suggests its incapability of utilizing extracellular polyamines. Nitrate and nitrite transporters have been found in the genome of strain E4A9. It processes genes involved in nitrate and nitrite reduction (*nasAB* and *nirBD*, respectively) and is lack of genes involved in denitrification, nitrogen fixation and anammox. Thus, nitrate and nitrite could act as electron acceptors to generate ammonium, subsequently being utilized by strain E4A9 as a reduced nitrogen source. The genome of *C. marinus* E4A9<sup>T</sup> is lack of urease (*ureABC*); however it harbors genes involved in urea decomposition, including urea carboxylase-related ABC transporter, urea carboxylase-related aminomethyltransferase, urea carboxylase and allophanate hydrolase, suggesting its capability of utilizing urea as a C or N source in the

**Table 5** Number of genes associated with general COG functional categories

Code	Value	%age <sup>a</sup>	Description
J	156	4.73	Translation, ribosomal structure and biogenesis
A	–	–	RNA processing and modification
K	190	5.76	Transcription
L	212	6.43	Replication, recombination and repair
B	1	0.03	Chromatin structure and dynamics
D	30	0.91	Cell cycle control, Cell division, chromosome partitioning
V	46	1.40	Defense mechanisms
T	168	5.10	Signal transduction mechanisms
M	193	5.86	Cell wall/membrane biogenesis
N	44	1.33	Cell motility
U	101	3.06	Intracellular trafficking and secretion
O	124	3.76	Posttranslational modification, protein turnover, chaperones
C	228	6.92	Energy production and conversion
G	187	5.67	Carbohydrate transport and metabolism
E	220	6.67	Amino acid transport and metabolism
F	64	1.94	Nucleotide transport and metabolism
H	146	4.43	Coenzyme transport and metabolism
I	199	6.04	Lipid transport and metabolism
P	174	5.28	Inorganic ion transport and metabolism
Q	111	3.37	Secondary metabolites biosynthesis, transport and catabolism
R	413	12.53	General function prediction only
S	289	8.77	Function unknown
–	770	23.36	Not in COGs

<sup>a</sup>The total is based on the total number of protein coding genes in the genome



environment [27]. (iii) Sulfur metabolism. Strain E4A9<sup>T</sup> possesses genes involved in assimilatory sulfate reduction (e.g. *cysND*, *cysC*, *cysH*, *cysJI*). Sulfate can be reduced to sulfide, subsequently being incorporated into amino acids. Genes involving in alkanesulfonate assimilation (arylsulfatase and FMN reductase) are present in the genome of strain E4A9, suggesting its capability of utilizing organic sulfur compounds. However, it missed transporter genes for the uptake of extracellular alkanesulfonates. (iv) Phosphorus metabolism. Strain E4A9 is lack of genes for inorganic P storage as polyphosphate (*ppk*), as well as transport (*phnCDE*) and cleavage (*phnGHIJKN*) of organic P in the form of phosphonates [28]. While strain E4A9 possesses the high-affinity phosphate transport system (*pstSCAB*) and regulatory genes (*phoUBR*), indicating an alternative strategy for maintaining a reliable supply of phosphorus [29].

## Conclusions

The complete genome sequence of *C. marinus* E4A9<sup>T</sup> contains a circular chromosome as well as two large

circular plasmids and provides an insight into the genomic basis of its esterases production ability. Our data implies *C. marinus* E4A9<sup>T</sup> is a potential candidate in biotechnological application and facilitates the understanding for further industrial and biotechnological applications of esterases.

## Abbreviations

CDS: Coding sequence; CRISPRs: Clustered regularly interspaced short palindromic repeats; KAAS: KEGG automatic annotation server; KEGG: Kyoto encyclopedia of genes and genomes; ORF: Open reading frame

## Acknowledgments

This work was supported by grants from the National Natural Science Foundation of China (No. 41406174 and 31770004), the National Key Basic Research Program of China (2014CB441503) and the Natural Science Foundation of Zhejiang Province (LR17D060001).

## Authors' contributions

XX and CW organized the study. YW and YH performed laboratory experiments. YW, HC and LX analyzed the data. YW drafted the manuscript. XX and QL edited the manuscript. All authors read and approved the final manuscript.

**Competing interests**

The authors declare that they have no competing interests.

**Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 15 October 2017 Accepted: 5 December 2017

Published online: 21 December 2017

**References**

- Jiang X, Xu X, Huo Y, Wu Y, Zhu X, Zhang X, et al. Identification and characterization of novel esterases from a deep-sea sediment metagenome. *Arch Microbiol.* 2012;194:207–14.
- Lopez-Lopez O, Cerdan ME, Gonzalez Siso MI. New extremophilic lipases and esterases from metagenomics. *Curr Protein Pept Sci.* 2014;15:445–55.
- Argigny JL, Jaeger KE. Bacterial lipolytic enzymes: classification and properties. *Biochem J.* 1999;343:177–83.
- Bornscheuer UT. Microbial carboxyl esterases: classification, properties and application in biocatalysis. *FEMS Microbiol Rev.* 2002;26:73–81.
- Yang S, Qin Z, Duan X, Yan Q, Jiang Z. Structural insights into the substrate specificity of two esterases from the thermophilic *Rhizomucor miehei*. *J Lipid Res.* 2015;56:1616–24.
- Xu XW, Wu YH, Wang CS, Wang XG, Oren A, Wu M. *Croceicoccus marinus* gen. Nov., sp. nov., a yellow-pigmented bacterium from deep-sea sediment, and emended description of the family *Erythrobacteraceae*. *Int J Syst Evol Microbiol.* 2009;59:2247–53.
- Lee K-B, Liu C-T, Anzai Y, Kim H, Aono T, Oyaizu H. The hierarchical system of the 'Alphaproteobacteria': description of *Hyphomonadaceae* fam. Nov., *Xanthobacteraceae* fam. Nov. and *Erythrobacteraceae* fam. Nov. *Int J Syst Evol Microbiol.* 2005;55:1907–19.
- Huang Y, Zeng Y, Feng H, Wu Y, Xu X. *Croceicoccus naphthovorans* sp. nov., a polycyclic aromatic hydrocarbons-degrading and acylhomoserine-lactone-producing bacterium isolated from marine biofilm, and emended description of the genus *Croceicoccus*. *Int J Syst Evol Microbiol.* 2015;65:1531–6.
- Wu YH, Li GY, Jian SL, Cheng H, Huo YY, Wang CS, et al. *Croceicoccus pelagius* sp. nov. and *Croceicoccus mobilis* sp. nov., isolated from marine environments. *Int J Syst Evol Microbiol.* 2016;66:4506–11.
- Yabuuchi E, Kosako Y. Order IV. *Sphingomonadales* ord. Nov. In: Garrity GM, Brenner DJ, Krieg NR, Staley JT, editors. *Bergey's manual of systematic bacteriology*, vol. Volume 2, Part C. Second ed. New York: Springer; 2005. p. 230–3.
- Euzéby J. Validation list no. 107. List of new names and new combinations previously effectively, but not validly, published. *Int J Syst Evol Microbiol.* 2006;56:1–6.
- Garrity GM, Bell JA, Lilburn T. Class I. *Alphaproteobacteria* class. Nov. In: Garrity GM, Brenner DJ, Krieg NR, Staley JT, editors. *Bergey's manual of systematic bacteriology*, Second Edition, Volume 2, Part C. New York: Springer; 2005. p. 1.
- Garrity GM, Bell JA, Lilburn T. Phylum XIV. *Proteobacteria* phyl. Nov. In: Garrity GM, Brenner DJ, Krieg NR, Staley JT, editors. *Bergey's manual of systematic bacteriology*, Second Edition, Volume 2, Part B. New York: Springer; 2005. p. 1.
- Lagesen K, Hallin P, Rodland EA, Staerfeldt HH, Rognes T, Ussery DW. RNAMmer: consistent and rapid annotation of ribosomal RNA genes. *Nucleic Acids Res.* 2007;35:3100–8.
- Lowe TM, Chan PP. tRNAscan-SE on-line: integrating search and context for analysis of transfer RNA genes. *Nucleic Acids Res.* 2016;44:W54–W7.
- Aziz RK, Bartels D, Best AA, DeJongh M, Disz T, Edwards RA, et al. The RAST server: rapid annotations using subsystems technology. *BMC Genomics.* 2008;9:75.
- Tatusov RL, Natale DA, Garkavtsev IV, Tatusova TA, Shankavaram UT, Rao BS, et al. The COG database: new developments in phylogenetic classification of proteins from complete genomes. *Nucleic Acids Res.* 2001;29:22–8.
- Finn RD, Coghill P, Eberhardt RY, Eddy SR, Mistry J, Mitchell AL, et al. The Pfam protein families database: towards a more sustainable future. *Nucleic Acids Res.* 2016;44:D279–85.
- Petersen TN, Brunak S, von Heijne G, Nielsen H. SignalP 4.0: discriminating signal peptides from transmembrane regions. *Nat Methods.* 2011;8:785–6.
- Krogh A, Larsson B, von Heijne G, Sonnhammer EL. Predicting transmembrane protein topology with a hidden Markov model: application to complete genomes. *J Mol Biol.* 2001;305:567–80.
- Grissa I, Vergnaud G, Pourcel C. CRISPRFinder: a web tool to identify clustered regularly interspaced short palindromic repeats. *Nucleic Acids Res.* 2007;35:W52–7.
- Kanehisa M, Araki M, Goto S, Hattori M, Hirakawa M, Itoh M, et al. KEGG for linking genomes to life and the environment. *Nucleic Acids Res.* 2008;36:D480–4.
- Moriya Y, Itoh M, Okuda S, Yoshizawa AC, Kanehisa M. KAAS: an automatic genome annotation and pathway reconstruction server. *Nucleic Acids Res.* 2007;35:W182–5.
- Stothard P, Wishart DS. Circular genome visualization and exploration using CGView. *Bioinformatics.* 2005;21:537–9.
- Castellani A, Chalmers AJ. *Manual of tropical medicine*. 3rd ed. New York: Williams Wood and Co; 1919.
- Skerman VDB, McGowan V, Sneath PHA. Approved lists of bacterial names. *Int J Syst Evol Microbiol.* 1980;30:225–420.
- Kanamori T, Kanou N, Atomi H, Imanaka T. Enzymatic characterization of a prokaryotic urea carboxylase. *J Bacteriol.* 2004;186:2532–9.
- Moran MA, Belas R, Schell MA, Gonzalez JM, Sun F, Sun S, et al. Ecological genomics of marine Roseobacters. *Appl Environ Microbiol.* 2007;73:4559–69.
- Liu Q, Wu YH, Cheng H, Xu L, Wang CS, Xu XW. Complete genome sequence of bacteriochlorophyll-synthesizing bacterium *Porphyrabacter neustonensis* DSM 9434. *Stand Genomic Sci.* 2017;12:32.
- Field D, Garrity G, Gray T, Morrison N, Selengut J, Sterk P, et al. The minimum information about a genome sequence (MIGS) specification. *Nat Biotechnol.* 2008;26:541–7.
- Woese CR, Kandler O, Wheelis ML. Towards a natural system of organisms: proposal for the domains Archaea, bacteria, and Eucarya. *Proc Natl Acad Sci U S A.* 1990;87:4576–9.
- The Gene Ontology C, Ashburner M, Ball CA, Blake JA, Botstein D, Butler H, et al. Gene ontology: tool for the unification of biology. *Nat Genet.* 2000;25:25–9.

Submit your next manuscript to BioMed Central and we will help you at every step:

- We accept pre-submission inquiries
- Our selector tool helps you to find the most relevant journal
- We provide round the clock customer support
- Convenient online submission
- Thorough peer review
- Inclusion in PubMed and all major indexing services
- Maximum visibility for your research

Submit your manuscript at  
www.biomedcentral.com/submit

